




WE CARE LIKE FAMILY




Air/Water Syringe Tips

Do not require an adapter and can be changed easily

- Fit most air/water syringes; no adapter required
- Eliminate air and water contamination
- Superior infection control
- Angled tip provides ideal placement of the tip anywhere in the patient's mouth
- These products are disposable
- Should not be autoclaved



See our website for our Pre-Bent Flow Tips (#2211-2214)



Reorder No.	Description	Quantity
2231	Blue/Clear	40/250/Cs
2232	Clear/Clear	40/250/Cs
2233	White/Clear	40/250/Cs



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

© 2016 Dynarex Corporation