


Standard Biopsy Punch

Integra® Miltex® offers an extensive line of sterile disposable biopsy punches. Each disposable biopsy punch features a seamless stainless steel design and a ribbed handle to help provide comfort and control. Sizes are clearly marked on each punch.

- Sterile
- Ten sizes for a wide range of applications
- Assortment box contains 10 each of the following sizes: 2 mm, 3 mm, 4 mm, 5 mm, and 6 mm

Reference	Size	Quantity
33-31AA	1.0 mm	50/box
33-31A	1.5 mm	50/box
33-31	2.0 mm	50/box
33-31B	2.5 mm	50/box
33-32	3.0 mm	50/box
33-33	3.5 mm	50/box
33-34	4.0 mm	50/box
33-35	5.0 mm	50/box
33-36	6.0 mm	50/box
33-37	8.0 mm	50/box
33-38	Assorted sizes	50/box




Surgical Blades

Carbon, Stainless Steel or Silicone Coated

CS* Reference	SS** Reference	SC*** Reference	Size	Quantity
4-110	4-310	4610	10	100/box
4-111	4-311	4611	11	100/box
4-112	4-312	-	12	100/box
4-112B ^Δ	4-312B ^Δ	-	12B ^Δ	100/box
4-115	4-315	4615	15	100/box
4-115C ^{ΔΔ}	4-315C ^{ΔΔ}	-	15C ^{ΔΔ}	100/box
4-120	4-320	-	20	100/box
4-121	4-321	-	21	100/box
4-122	4-322	-	22	100/box
4-123	4-323	-	23	100/box
4-125	4-325	-	25	100/box

- Sterile
- Razor-sharp cutting edge
- Sizes clearly marked on each blade
- Individually packaged

*Carbon Steel ** Stainless Steel *** Silicone Coated
^ΔB is a double edged blade ^{ΔΔ}C is an angled blade



Surgical Blade Remover

Reference	Description
4-100	Blade Safe Surgical Blade Remover For Any Style Handle



4-100

Standard Scalpels

- Sterile
- Stainless steel, razor-sharp cutting edge
- Safety sleeve for blade
- Individually packaged

Reference	Size	Quantity
4-410	10	10/box
4-411	11	10/box
4-412	12	10/box
4-415	15	10/box
4-420	20	10/box
4-421	21	10/box
4-422	22	10/box
4-423	23	10/box



4-410



4-411



4-412



4-415



4-420



4-421



4-422




4-423


Safety Scalpels

- Sterile
- Stainless steel, razor-sharp cutting edge
- Metric ruler on each scalpel
- Individually packaged

Reference	Size	Quantity
4-510	10	10/box
4-511	11	10/box
4-512	12	10/box
4-515	15	10/box
4-515C	15C Angled Blade	10/box
4-522	22	10/box



4-510




4-511




4-512



4-515



4-515C



4-522

Directions For Use

Gently push slider to right and completely forward in one simple motion. Blade locks into place with audible confirmation. To retract, push slider to right and back to return blade to original protected position.

