

ADJUSTABLE GEL-PACK HEAD POSITIONER

Reorder #	Description	Size	Wt.	Unit
<input type="checkbox"/> 914364	Adjustable Positioner	17"x7"x7"	1.4 lbs	Ea.
<input type="checkbox"/> 503099 medium	Replacement gel pk (2 required)	4"x 8.5"	.5 lbs	Ea.

PURPOSE

Individuals who require assistance in upper body support typically have difficulty in maintaining symmetrical head positioning. The one size fits all head rests often are uncomfortable and create an environment for ear irritation and consequently, sores on the ears. This adjustable head support accommodates a wide range of head sizes and with the utilization of a cooling gel will prevent most cases of sores on the ears. This universal head support controls lateral head movement and with adjustable straps, fits any high back chair.

USE INSTRUCTIONS

Utilizing the adjustable straps centrally locate and affix the headrest to the chair (illustration 1). Raise or lower the head rest to allow for maximum comfort and torso/head positioning. One side of the foam wedge is removable and adjustable to maximize comfort (illustration 2). The unit is provided with two pockets to house the gel for ear protection (illustration 3) which can be removed if desired.

CARE INSTRUCTIONS

The adjustable head rest is covered with wipe clean vinyl and can be cleaned with a mild cleaner. Do not launder this product.

