



Medication preparation |

Medication delivery |



Sharps disposal

Cardinal Health™ Monoject™ Insulin and Tuberculin Safety Syringes

Safety needles designed to protect

for strong performance and value



Protect your staff on their journey to the patient bedside

When you standardize with Cardinal Health™ Sharps Safety Products, you get a line of safety-engineered products designed to decrease the risk of needlestick injuries and minimize exposure to blood and bodily fluid. In fact, the CDC estimates that up to 75% of sharps injuries can be prevented simply by choosing safer medical devices.¹

Choose Monoject™ Insulin and Tuberculin Safety Syringes

Monoject™ Insulin and Tuberculin Safety Syringes utilize a barrel-based safety technology designed to keep patients and clinicians safe with increased flexibility.

The exclusive design allows staff to travel safely to a patient's bedside without the risk of an exposed needle.


- 1 Full coverage safety shield, not an accessory
- 2 True "transport mode"
- 3 Two-click activation
- 4 Clear barrel and bold, easy-to-read graduations



Monoject™ Insulin and Tuberculin Safety Syringes

Monoject™ Insulin Safety Syringes with Permanent Needles


- Sterile, single use and individually packaged
- Soft pack packaging
- Low dead space
- Permanently attached needle
- Flat plunger tip
- For U-100 insulin
- Supports diabetes care
- Not manufactured with DEHP or natural rubber latex



Order code	Description	Syringe volume	Needle gauge	Needle length	Box qty.	Case qty.
8881511144	⅓ mL 29 G x ½ in. (0.33 x 13 mm)	⅓ mL	29 G	½ in.	100	500
8881511136	½ mL 29 G x ½ in. (0.33 x 13 mm)	½ mL	29 G	½ in.	100	500
8881511110	1 mL 29 G x ½ in. (0.33 x 13 mm)	1 mL	29 G	½ in.	100	500
8881511344	⅓ mL 30 G x ⅝ in. (0.3 x 8 mm)	⅓ mL	30 G	⅝ in.	100	500
8881511336	½ mL 30 G x ⅝ in. (0.3 x 8 mm)	½ mL	30 G	⅝ in.	100	500
8881511310	1 mL 30 G x ⅝ in. (0.3 x 8 mm)	1 mL	30 G	⅝ in.	100	500

Monoject™ Tuberculin Safety Syringes with Permanent Needles

- Sterile, single use and individually packaged
- Soft pack packaging
- Low dead space
- Permanently attached needle
- Flat plunger tip
- Graduations: 0.01 mL
- Not manufactured with DEHP or natural rubber latex



Order code	Description	Syringe volume	Needle gauge	Needle length	Box qty.	Case qty.
8881511235	1 mL 25 G x ⅝ in. (0.5 x 16 mm)	1 mL	25 G	⅝ in.	100	500
8881511201	1 mL 28 G x ½ in. (0.36 x 13 mm)	1 mL	28 G	½ in.	100	500

For more information or to order Monoject™ Sharps Safety Products, contact your sales representative, call **800.964.5227** or visit **cardinalhealth.com/monoject**

Access the Monoject™ Insulin and Tuberculin Safety Syringes instructions for use at **cardinalhealth.com/monoject/IFU**

References: 1. NIOSH. Alert: preventing needlestick injuries in health care settings. Cincinnati, OH: US Dept of Health and Human Services; 1999. <https://www.cdc.gov/niosh/docs/2000-108/pdfs/2000-108.pdf>.

© 2021 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO and MONOJECT are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other trademarks are the property of their respective owners. Patent cardinalhealth.com/patents. Lit. No. 2GM21-1652238-02 (09/2021)

